

Do You Know What You're Made Of?

Go beyond the scale and get an InBody Total Body Composition Assessment from WPV Inc.!





THE TEST

An InBody Test is easy, quick, and precise. Simply stand on the device and grab the hand electrodes to begin.

THE METHOD

The hand and foot electrodes send safe currents through your body and measure the impedance values.

THE PROCESS

From these measured impedance values, your body composition is determined without the use of statistical data.

THE RESULTS

In 15 seconds, the InBody will finish analyzing your entire body and print out your full-page results sheet.

Learn more by visiting: www.WPVinc.com

